Art. 85153

85153 05-2021 Ver. 2.0

**Gebrauchsanweisung**  
**Instruction manual**  
**Gebruiksaanwijzing**  
**Mode d'emploi**  
**Manual de instrucciones**  
**Istruzioni per l'uso**  
**Manual de instruções**  
**Εγχειρίδιο οδηγιών χρήσης**  
**Käyttöohje**  
**Bruksanvisning**  
**Bruksanvisning**  
**Brugsanvisning**  
**Instrukcja obsługi**  
**Návod k použití**  
**Návod na použitie**  
**Navodila za uporabo**  
**Uputa o uporabi**  
**Használati útmutató**  
**Instrucțiuni de utilizare**  
**Инструкция за употреба**  
**Kasutusjuhend**  
**Lietošanas instrukcija**  
**Naudojimo instrukcija**  
**Руководство по использованию**

medisana GmbH  
 Carl-Schurz-Str. 2  
 41460 NEUSS  
 GERMANY

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**EN**  
**IMPORTANT INFORMATION!**  
**BE SURE TO RETAIN FOR YOUR REFERENCE!**

- 1 On/off switch (Clockwise/anti-clockwise rotation)
- 2 Operating Display
- 3 „-“ Button to decrease rotation speed
- 4 „+“ Button to increase rotation speed
- 5 Holder for attachments
- 6-12 attachments (in the order depicted: sapphire cone, felt cone, sapphire disc (rough), sapphire disc (fine), cylindrical burr, flame burr, callus remover)

**Safety instructions**

- Never cover the device.
- Never allow the device to operate unattended.
- The device is only fully disconnected from the power supply when the power adaptor has been removed from the power socket.
- Always remove the power adaptor when you are no longer using the device. This is particularly important even when just taking a break.
- Before connecting the device to your power supply, make sure that the supply voltage specified on the data plate is compatible with your mains supply.
- Use the device only for its intended purpose and as directed within this instruction manual.
- The device is not intended for commercial use.
- Do not use any auxiliary parts that have not been recommended by the manufacturer.
- Avoid contact with pointed or sharp objects. The warranty shall be rendered null and void if the device is not used for its intended purpose.
- This device may be used by children over 8 years of age as well as persons with reduced physical, sensory or mental capacities or lack of experience and knowledge, providing they are supervised or have been instructed regarding the safe use of this device and understand potential hazards involved with said use.
- Children must be supervised to ensure they do not play with the device.
- Cleaning and user maintenance must not be performed by children unless they are supervised.
- Do not attempt to repair the device yourself in the event of any defects. This will not only void any warranty claim, but may also create a serious hazard (fire, electric shock, injury). Only have repairs carried out at authorised service points.
- Only plug in the power adaptor when the device is switched off.
- Do not reach for the device if it has fallen into water. Remove the power adaptor from the socket immediately.
- Keep the power supply cable away from hot surfaces and open flames.
- Do not use the device in rooms in which aerosols (sprays) are being used or where pure oxygen is being administered.
- Make sure to arrange the mains cable in a way that does not pose a trip or asphyxiation hazard.

- Never use any electrical device while bathing or showering. Do not store or keep any electrical devices in places where they might fall into a bath or a sink.
- Do not use the device if the power cable has been damaged, the power adaptor is defective, the device does not work perfectly or if it has fallen into water or been damaged.
- Never carry, pull on or rotate the device by the mains cable and ensure the cable is not subject to entrapment.
- Avoid overheating due to excessive use.
- Uninterrupted treatment should not take any longer than 15 minutes. If this time is exceeded, switch the device off for at least half an hour and allow it to cool down.
- Only pull the power adaptor out of the socket using dry hands, when the device is switched off and no longer in use as well as when you need to replace or clean the attachments.

**Important recommendations for your health**

- If diagnosed with diabetes or other illnesses, contact your family doctor before attempting to use this device.
- Pregnant women should follow necessary precautionary measures and consider their own physical state. Consult a doctor, if necessary, before use.
- Do not use the device on any part of the body exhibiting swelling, burns, inflammations, skin rashes, wounds or other sensitive areas.
- The treatment should feel pleasant. Discontinue the treatment and consult your doctor if you feel any pain or find the treatment unpleasant.
- Do not soak the skin of your hands and feet before use.

**How to use**

**WARNING**  
**Risk of injury! Excessively long or intensive application can damage deeper skin tissue.**

**On/off switch**

To switch the device on, push the switch 1 in one of two given directions. Push the switch downwards (for clockwise rotation direction, operating display 2 will light up green) or upwards (for anti-clockwise rotation direction, operating display 2 will light up red).

To switch the device off, push the switch back towards the middle (O).

**+ / - Buttons to increase/decrease rotation speed**

The **MP 815** features two different rotation speeds (level 1 and level 2). The device starts at a low rotation speed once switched on (level 1). Press 4 to increase the rotation speed (level 2). Press 3 to decrease the rotation speed (level 1). Each time it is used, the device should be started at a low rotation speed, which you can increase slowly to meet your individual requirements.

**Changing the attachments**

The device is equipped with a hexagonal attachment holder. First, turn the device off. Do not just pull out the attachment you do not need. Place the new attachment into the holder and turn it slightly to ensure it clicks into position.

**Attachments**

Designated use of all attachments (in the order depicted):

- **Sapphire cone** 6: To remove dry hard skin.
- **Felt cone** 7: To smoothen the edge of the nail and clean the nail surface.
- **Sapphire disc, rough** 8, **fine** 9: To shorten the nails by filing in the direction of growth. Hold the device and run the nail along the disc. File in one direction only!
- **Cylindrical burr** 10: To file down hardened nail surfaces.
- **Flame burr** 11: To carefully remove ingrown nails.
- **Callus remover** 12: To remove hard skin or calluses on the sole or heel of the foot.

**Cleaning and Care**

- The power adaptor for all electrical devices should always be removed from the socket.
- Never use aggressive cleaning agents or stiff brushes.
- To clean the device, use a soft cloth that has been dampened with a mild soap solution or water. Never immerse the device in water when cleaning and ensure that water does not otherwise enter the device. Only use the device again when it is completely dry.
- The attachments can be cleaned with a brush (e.g. an old toothbrush) for hygiene reasons. It is also recommended that the grinding surfaces are cleaned with commercially available disinfectant spray or 90% alcohol to prevent the transmission of infection. The attachment pins must be kept free from grease. Let the attachments dry completely before (re-)attaching them to the device.
- Store the device in a cool, dry place.
- Ensure the power cable is neatly packed to avoid cable breaks.

**Instructions for disposal**

Do not throw this device away with ordinary household waste. Consumers are obliged to bring all electric or electronic devices, irrespective of whether they contain harmful substances or not, to a community or retailer's collection site so they can be transferred to the proper channels for disposal. Contact your local authority or retailer regarding disposal.

**Technical data**

Name and model: medisana Manicure/ Pedicure Device MP 815; Rotation speed: approx. 3.800 U/Min. (Level 1); approx. 4.600 U/Min. (Level 2); Dimensions: approx. 15x3.8x3.7cm; Weight: approx. 942g; Item number: 85153; EAN number: 4015588851537; Electrical data: Comprehensive data regarding the adaptor can be taken from the data sheet provided.

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